

Optimizing Outcomes with Standardized Tracheostomy and Laryngectomy HME Use

Carmin Bartow, MS, CCC-SLP | Casey Richardson, MA, CCC-SLP
Senior Clinical Educators – Atos Medical

Why Pulmonary Health Matters

Breathing through a tracheostomy or laryngectomy results in inhaled air that is:

- Dry
- Cold
- Unfiltered

Consequences may include:

- Thickened secretions
- Increased coughing
- Airway inflammation and narrowing
- Reduced ventilation and pulmonary health

Artificial Humidification Options

Active External Humidification (delivered via tubing and a tracheostomy mask)

Disadvantages include:

- Increased noise and sleep disruption
- Reduced mobility
- Requires electricity and ongoing maintenance
- Higher cost and training burden
- Lower patient satisfaction and adherence

Passive Humidification: heat and moisture exchangers (HMEs)

- Capture heat and moisture from exhaled air
- Release warmth and humidity during inhalation
- Provide a barrier to gross airborne particles

Evidence Supporting HME Use: Laryngectomy

Level 1 evidence supports HMEs over external humidification post-laryngectomy
Consistent HME use (24/7):

- Reduces coughing and forced expectoration
- Improves secretion management
- Decreases mucus plugging and in-hospital complications

RCTs and cohort studies demonstrate:

- Fewer pulmonary symptoms
- Reduced suctioning needs
- Improved comfort and mobility
- Strong nurse and patient preference

Provox Life HMEs:

- Further reduce pulmonary symptom burden
- Improve quality-of-life measures (CASA-Q)
- Optimize performance in everyday environments

Evidence Supporting HME Use: Tracheostomy

HMEs are effective and well tolerated in spontaneously breathing patients

Hygroscopic foam HMEs (e.g., calcium chloride-treated):

- Provide superior moisture retention
- Improve secretion consistency

TrachPhone HME outcomes:

- 97% immediate post-operative tolerance
- Improved communication
- Reduced noise, suctioning, and training needs
- Cost savings and improved workflow efficiency

Driving Practice Change

Evidence supports standardizing HME use across care settings

Successful implementation includes:

- Interdisciplinary engagement
- Education and workflow alignment
- Barrier identification and removal
- Ongoing reinforcement and sustainability planning

Implementation Support

Laryngectomy Pulmonary Kits & Care Pathways

- Promote consistent, long-term HME adherence
- Reinforce best practices across settings

TrachPhone Implementation Program

- Guided transition from external humidification to heat and moisture exchangers (HMEs)
- Reduces practice variability
- Simplifies education for staff, patients, and caregivers
- Integrates evidence into daily clinical workflows

Key Takeaways

- Surgical airways significantly impair natural humidification and filtration
- HMEs are an evidence-based standard of care for pulmonary health
- Consistent HME use improves symptoms, reduces complications, and enhances quality of life
- Standardization supports better outcomes and more efficient care delivery

References

- Bień S, Okła S, van As-Brooks CJ, Ackerstaff AH. The effect of a heat and moisture exchanger (Provox HME) on pulmonary protection after total laryngectomy: a randomized controlled study. *Eur Arch Otorhinolaryngol*. 2010;267(3):429–435. doi:10.1007/s00405-009-1018-4.
- de Kleijn BJ, van As-Brooks CJ, Wedman J, van der Laan BFAM. Clinical feasibility study of ProTrach DualCare, a new speaking valve with heat and moisture exchanger for tracheotomized patients. *Laryngoscope Invest Otolaryngol*. 2017;2(6):453–458.
- Foreman A, De Santis RJ, Sultanov F, Enepekides DJ, Higgins KM. Heat and moisture exchanger use reduces in-hospital complications following total laryngectomy: a case-control study. *J Otolaryngol Head Neck Surg*. 2016;45(1):40. doi:10.1186/s40463-016-0154-2.
- Intensive Care Society Standards. Tracheostomy Care. London (UK): Intensive Care Society; 2014. Available from: <https://ics.ac.uk/asset/E0B8C2B2-2B70-4A18-839EB859969DBB74/> Kearney A, Norris K, Bertelsen C, Samad I, Cambridge M, Croft G, et al. Adoption and utilization of heat and moisture exchangers (HMEs) in the tracheostomy patient. *Otolaryngol Head Neck Surg*. 2023.
- Kotter JP. *Leading change*. Boston (MA): Harvard Business Review Press; 2012.
- Merol JC, Charpiot A, Langagne T, Hemar P, Ackerstaff AH, Hilgers FJ. Randomized controlled trial on postoperative pulmonary humidification after total laryngectomy: external humidifier versus heat and moisture exchanger. *Laryngoscope*. 2012;122(2):275–281.
- National Tracheostomy Safety Project. Day-to-day management of tracheostomies and laryngectomies: humidification. 2014. Available from: <https://tracheostomy.org.uk/storage/files/HumidificationNew.pdf>
- Parrilla C, Minni A, Bogaardt H, Macri GF, Battista M, Roukos R, et al. Pulmonary rehabilitation after total laryngectomy: a multicenter time-series clinical trial evaluating the Provox XtraHME in HME-naïve patients. *Ann Otol Rhinol Laryngol*. 2015;124(9):706–713. doi:10.1177/0003489415579219.
- van den Boer C, Vas Nunes JH, Muller SH, van den Brekel MW, Hilgers FJ. Water uptake performance of hygroscopic heat and moisture exchangers after 24-hour tracheostoma application. *Otolaryngol Head Neck Surg*. 2014;150(6):999–1004.
- Vitacca M, Clini E, Foglio K, et al. Hygroscopic condenser humidifiers in chronically tracheostomized patients who breathe spontaneously. *Eur Respir J*. 1994;7(11):2026–2032.
- Ward EC, Hancock K, Boxall J, Burns CL, Spurgin AL, Lehn B, et al. Post-laryngectomy pulmonary and related symptom changes following adoption of an optimal day-and-night heat and moisture exchanger regimen. *Head Neck*. 2023;45(4):939–951. doi:10.1002/hed.27323.

Atos Medical Inc

5000 S. Towne Drive, Suite 200, New Berlin, WI 53151
+1 800 217 0025 | info.us@atosmedical.com | www.atosmedical.us

Atos and the Atos Medical logo are trademarks of Coloplast A/S.
© 2026 Coloplast A/S. All rights reserved. | MC-3412, 021726