# Tracoe Purofoam-T

## Non-adhesive tracheostomy dressing

## Instructions for use

### Tracoe Purofoam-T

#### Non-adhesive tracheostomy dressing

Tracoe Purofoam-T can be used by lay persons and/or health care professionals. Read the instructions for use prior to using the product and retain for your information. If you are unsure how to use this product, or whether it is suitable for you, please consult a healthcare professional.

**Product Description:**  
A soft, sterile, conforming, non-adherent absorbent foam dressing with an aesthetically pleasing semi-permeable film backing and key-hole fenestration to allow ease of placing around a tube, cannula or peg site.

The dressing offers softness and air permeability, perfect for the treatment of low/medium (Tracoe Purofoam-T Large Lite) and medium/high (Tracoe Purofoam-T Large/Small) exuding wounds and secretions.

The dressing provides a cushioning effect surrounding the wound. It is absorbent, lightweight and provides a rapid uptake of exudate. The dressing has a bacteria-proof polyurethane backing, is breathable and prevents strike-through and maceration.

**Product Material Composition:**  
Polyurethane foam

**Intended Purpose/Indications for Use:**

Should be used for managing secretions around intact skin and improving comfort for patients with a tracheostomy tube, cannula or peg site.

**Intended Population/User:**

Tracoe Purofoam-T have no specific intended patient population, other than patients suffering a relevant medical indication.

**Contra-Indications:**

Do not use with oxidising agents such as hypochlorite solutions or hydrogen peroxide as these can break down the absorbent component of the dressing.

Not suitable for arterial bleeds, heavy bleeding wounds.

**Warnings and Precautions:**

Tracoe Purofoam-T is a single use product. Re-use of this dressing after being placed on the patient may result in patient infection and/or cross contamination. Re-sterilisation, reprocessing, cleaning and/or disinfection may also compromise the product characteristics, resulting in trauma or infection to the patient.

Opening the dressing pack compromises the sterile barrier. Any unused dressing should be discarded.

**Instructions for Use:**

Ensure the area has been cleaned. Gently dry the surrounding skin.

Apply with care – pink film side up and white side next to skin. Open the key-hole and gently place the dressing around the tube, rod or cannula. Training may be required to avoid potential trauma. If in doubt, seek professional medical advice from a doctor or healthcare professional.

Tracoe Purofoam-T can be left in place for up to 7 days but should be changed if soiled or damp from secretions or exudate. If exudate

Mild to Moderate:  
• Delayed wound healing  
• Further tissue damage  
• Infection  
• Leakage  
• Necrosis  
• Shock  
• Wound bed becoming dry  
**In severe cases:**  
• Anaphylactic shock  
• Sepsis  
• Scarring

**Device Supplied Sterile:**  
• Store devices in a manner that provides protection from dust, light, insects, vermin and extreme temperature and humidity.

• Device is supplied sterile. When provided sterile, always store the product in the original protective packaging.

**Handling and Storage:**  
• Store at room temperature.  
• Keep dry.  
• Keep out of direct sunlight.

**Disposal:**

Dispose of hygienically after use.  
Dispose as clinical waste in accordance with local regulations if used in an acute or hospital setting.

**Warranty:**

All warranty rights are lost if changes or modifications are carried out by an unauthorised personnel.

The manufacturer does not take responsibility for any effects on safety, reliability or performance of the product if the product is not used in conformity with the instructions for use.

**Notice to EU Users and/or Patients (Complaints/Incidents and Information):**

**Any serious incident that has occurred in relation to the device, should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.**

Any other grounds for dissatisfaction relating to the quality of the product should be notified to Brightwake Ltd or its representative.

For all complaints, please give the name and reference along with the batch number of the component(s), your name and address and an exhaustive description of the event to help Brightwake Ltd understand the cause of the complaint. When possible, please send the device to Brightwake Ltd for the investigation as per the following:  
complaints@brightwake.co.uk

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