PROVOX Dilator

Rx
ONLY

CE 2460

Provox Dilator 17/20

REF 8342, 8343


Figure 1


Figure 2


Figure 3


Figure 4


Figure 5


Figure 6


Figure 7


Figure 8


Figure 9


Figure 10


Figure 11

Figure 12


Figure 13


Figure 14


Figure 15

Prescription Information

CAUTION: United States Federal law restricts this device to sale, distribution, and use by or on order of a physician or a licensed practitioner. The availability of this product without prescription outside the United States may vary from country to country.

Disclaimer

Atos Medical offers no warranty—neither expressed nor implied—to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

Contents

EN - ENGLISH	6
DE - DEUTSCH	13
FR - FRANÇAIS	22
IT - ITALIANO	30
ES - ESPAÑOL	38
PT - PORTUGUÊS	46
ET - EESTI	54
HU - MAGYAR	61
PL - POLSKI	68
HR - HRVATSK	76
EL - ΕΛΛΗΝΙΚΑ	83
RU - РУССКИЙ	92
JA - 日本語版	101
KO - 한국어	109
ZHTW - 繁體中文	116
HE - עברית	128
AR - العربية	133

ENGLISH

1. Intended Use

Provox Dilator 17 and Provox Dilator 20 are intended to dilate or stent the tracheoesophageal (TE) puncture in laryngectomized patients.

Intended User Group

The device is intended to be used by healthcare professionals, adult patients (18+ years), and lay caregivers.

2. Contraindications

- The device is not intended to be used by patients and/or lay caregivers who have not received adequate training from their clinician and, as assessed by a clinician who prescribes the device, have demonstrated the ability to understand and consistently follow the Instructions for Use (IFU) without clinician supervision.
- The device is not intended to be used for TE puncture dilation at the time of surgical creation of the TE puncture.

3. Description of the Device

Provox Dilator 17/20 is for single patient, multiple use only. The dilators are tapered curved silicone rods used for dilating (increase the diameter of) or stenting (blocking) TE punctures. Provox Dilator 17 shall be used with voice prostheses with an outer diameter of 17 Fr and Provox Dilator 20 with an outer diameter of 20 Fr.

Provox Dilator 17/20 shall only be prescribed for patient use by clinicians trained in the care and rehabilitation of laryngectomized patients.

Provox Dilator 17/20 has a small retention collar to prevent the dilator from gliding back to the thinner section. It also has a safety strap with medallion that can be taped to the neck to keep the dilator in place and that is designed to reduce the risk of accidental aspiration of the device.

For dilating the TE puncture, the tip of the dilator is inserted in the open TE puncture and slowly progressed forward and down into the esophagus (food pipe) until the desired diameter of the TE puncture is reached. The gradual increase in the diameter of the device gradually dilates the TE puncture to facilitate the subsequent insertion of the voice prosthesis.

For stenting the TE puncture, the tip of the dilator is inserted in the open TE puncture and slowly progressed forward and down into the esophagus (food pipe) until it reaches the desired diameter for blocking the TE puncture. The presence of the dilator keeps the TE puncture open and prevents esophageal contents from entering the trachea (windpipe) and lungs.

Provox Dilator 17/20 (see Figure 1)

Medallion (1)
Safety strap (2)
Device handle (3)
Notch (4)
Retention collar (5)
Tip (6)

The size is marked on the medallion.

Dimensions

Provox Dilator 17	The diameter is 11.5 Fr at the tip and increases to 19 Fr.
Provox Dilator 20	The diameter is 13 Fr at the tip and increases to 22 Fr.
Length between tip and notch	60 mm (2.36 in.)
Length between retention collar and notch	20 mm (0.79 in.)

Recommended Water-Soluble Lubricants

For use with the dilator

K-Y Jelly
Surgilube
Xylocaine (Lidocaine) Jelly*

**Local anaesthetic*

Caution: Lubricants or anaesthetizing gels that are not water-soluble shall not be used at any time.

4. Warnings

- Do not use the device if it shows any signs of structural damage. A damaged device may break and parts of it can enter the trachea or esophagus and have to be removed immediately.
- Do not use lubricating or anaesthetizing gel if you have allergies related to these substances.
- Never remove the medallion from the device. The medallion is of the utmost importance for safety. It helps to reduce the risk of the device accidentally falling into the trachea or the esophagus.
- Never use excessive force to insert or remove the dilator from the TE puncture. If you are unable to insert or

remove the dilator with gentle force, immediately consult your clinician.

- To avoid contamination, follow appropriate hygienic procedures.
- Always ensure mechanical integrity before each use.
- The device is intended for single patient use only. Do not let anyone else use your dilator or use a dilator from someone else. Use of one device by multiple patients may cause serious infection.

5. Precautions

General

- It is advisable to lubricate the dilator before the procedure. Use only water-soluble lubricants in an appropriate amount. Do not use any kind of oil or ointment during the procedure. Oil-based lubricants (e.g. Vaseline) may weaken, damage, or destroy the device or may enter your windpipe and lungs which can result in severe coughing and/or lipoid pneumonia.
- Make sure that the diameter of the dilator is adequate for the TE puncture. Over-dilatation may result in leakage around the voice prosthesis.
- Never use excessive force to insert or remove the dilator from the TE puncture, as this may cause bleeding or a false route.
- Do not use a contaminated dilator. A contaminated dilator may cause infection and inflammation.
- Do not pull the safety medallion during removal (see Step 7 in Section 6.2 Operating Instructions).
- Do not use the dilator for too long a period of time as this may damage the esophagus.

Professional Use

- Patients on anticoagulant therapy or with bleeding disorders should be carefully evaluated for the risk of hemorrhage prior to dilatation of the TE puncture.
- Carefully assess the suitability of the device for patients with small tracheostomas. The device takes up space in the tracheostoma and may affect the patient's ability to breathe.

Patient/Caregiver Use

To reduce the risk of infection,

- make sure that your hands are thoroughly clean and dry before using your hands in the area of your tracheostoma and before handling your voice prosthesis and/or accessories.
- clean and disinfect the device after each use according to instructions (see Section 7 Cleaning and Disinfecting).
- if the dilator looks dirty or has been air dried in an area with a risk of contamination, the dilator should be both cleaned and disinfected before use. A risk of contamination could

exist if the device has been in contact with a pet, dropped on the floor or other gross contamination, or if someone has a respiratory infection.

- during hospitalisation, it is important to both clean and disinfect the dilator after use but also before use since there is an increased risk of infection and contamination. In a hospital, it is better to use sterile water for cleaning and rinsing, rather than drinking water.

Consult your clinician if you

- bleed from the TE puncture tract. Bleeding may be caused by damage to the TE puncture tract, damage to the esophageal wall, or by bleeding from esophageal vessels. Minor bleeding during dilatation is normal. If a minor bleeding does not stop, you should consult your clinician.
- accidentally dislodge the dilator. The dilator may be used for stenting the TE puncture during a longer period of time. Accidental dislodgement during these circumstances may be followed by shrinkage of the TE puncture and subsequent closure.
- experience an allergic reaction against the lubricating or anaesthetizing gel. Read the IFU accompanying the used lubricant.

Side-Effects

Possible side-effects of using Provox Dilator 17/20 are coughing, gagging, and discomfort when the device is inserted, in situ, and/or removed from the TE puncture.

Minor bleeding may occur due to the dilatation of the tissues of the TE puncture.

Users with smaller sized tracheostomas may temporarily experience increased breathing resistance when the device is in situ as it occupies space within the tracheostoma.

Users may also experience irritation and/or discomfort from the safety strap and medallion.

Notify your clinician if any of these side-effects occur.

6. How to Use

6.1 Preparation

Always handle the dilator with clean hands or gloves and avoid touching parts that will enter the puncture.

Prior to each use, check the integrity of the safety strap by slightly pulling the medallion with one hand while holding the device handle with the other hand (see Figure 2). Discard the device if it shows any signs of damage.

It is advisable to lubricate the dilator with some water-soluble or anaesthetizing gel prior to use (see Figures 4–5 and see Section 3 Description of the Device—Recommended Water-Soluble Lubricants).

6.2 Operating Instructions

- Make sure that you have a magnifying mirror and sufficient light to illuminate the tracheostoma.
- Avoid swallowing once the voice prosthesis is removed and the puncture is open since this may lead to aspiration of saliva and/or coughing. Keeping your mouth open can help you avoid swallowing for short periods.
- Do not clean, remove, or otherwise manipulate the voice prosthesis (or Provox Dilator 17/20) immediately after a meal since the stimulation of the TE puncture tract immediately after eating may lead to gagging and/or vomiting.
- Make sure that you have completely understood and fully practiced all procedures under the supervision of your clinician before you perform any procedure on your own.

Patient/Caregiver Use

1. Rinse the dilator in drinking water before use (see Figure 3).
2. If required, lubricate the tip of the dilator that will enter the TE puncture:
 - a. Put a small amount of water-soluble lubricating gel on a clean surface, e.g. a gauze (see Figure 4).
 - b. Slide the tip of the dilator through the gel so that the tip is covered with a thin layer of gel (see Figure 5).
3. Remove the voice prosthesis from the TE puncture before inserting the dilator. Make sure the TE puncture tract is empty and/or open.
4. Probe the TE puncture with the tapered tip of the dilator with the curve of the dilator directed downwards, and slowly and gently insert the device until you reach the notch in front of the handle (see Figures 6–8).
5. Keep the dilator in this position for the desired amount of time. The time needed to sufficiently dilate the TE puncture tract to facilitate voice prosthesis insertion varies from person to person.
6. Carefully tape the safety medallion to the skin while you wait (see Figure 9).
7. Once the dilator has been in the TE puncture tract for the desired amount of time, gently remove the dilator by pulling the handle (see Figure 10). Do not pull the safety medallion.
8. Insert the voice prosthesis using the standard procedure.

Troubleshooting

There might be several possible causes as to why you are unable to insert the dilator, such as wrong angle/direction when inserting, shrinkage of puncture, or deformation of the TE puncture tract. If this should occur, never use force to insert the dilator. Consult your clinician immediately.

If you have inserted the dilator too far into the TE puncture tract, you might have used the product in the wrong way or might have an enlarged puncture tract. If this occurs, gently try to retrieve the device by pulling the medallion. You should never use force. If you cannot retrieve the device, consult your clinician immediately.

Professional Use

1. Rinse the dilator in drinking water before use (see Figure 3).
2. If required, lubricate the tip of the dilator that will enter the TE puncture:
 - a. Put a small amount of water-soluble lubricating gel on a clean surface, e.g. a gauze (see Figure 4).
 - b. Slide the tip of the dilator through the gel so that the tip is covered with a thin layer of gel (see Figure 5).
3. Remove the voice prosthesis from the TE puncture.
4. Probe the TE puncture with the tapered tip of the dilator with the curve of the dilator directed downwards and slowly and gently insert the device until you reach the stop in front of the handle (see Figures 6–8).
5. Keep the dilator in this position for the desired amount of time. The time needed to sufficiently dilate the TE puncture tract to facilitate voice prosthesis insertion varies from patient to patient.
6. Carefully tape the safety medallion to the skin while you wait (see Figure 9).
7. Once the dilator has been in the puncture tract for the desired amount of time, gently remove the dilator by pulling the handle (see Figure 10). Do not pull the safety medallion.
8. Insert a voice prosthesis according to normal procedure.

7. Cleaning and Disinfecting

Clean and disinfect the dilator after each use.

1. Rinse the device with lukewarm drinking water for 30 s.
2. In a clean bowl, mix two drops of dish soap in 250 ml (8.5 fl oz) warm drinking water, 50–60 °C (122–140 °F). Place the dilator in the mixture until it is fully submerged for 1 min (see Figure 11)

Caution: Be careful not to burn yourself with the warm cleaning water.

3. Take the dilator out of the mixture and rub the dilator for 15 s to remove any remaining mucus and debris (see Figure 12). If it is not possible to remove all mucus and debris, repeat from Step 1.
4. Rinse it in lukewarm drinking water for 30 s or until all soap has been removed (see Figure 13).

5. Disinfect the dilator by fully submerging it in room temperature 3% hydrogen peroxide for 60 min, or 70% ethanol for 10 min, or 70% isopropyl alcohol for 10 min (see Figure 14).
6. After cleaning and disinfecting the device, place the dilator on clean gauze and leave it to air dry (see Figure 15). When the dilator is dry, store it at room temperature in a clean closed container. Make sure that the stored device is protected from both heat and direct sunlight.

Caution: Do not clean or disinfect by any other method as this might cause product damage and injury.

8. How Supplied

The device is manufactured in a cleanroom and supplied non-sterile.

9. Device Lifetime and Disposal

Caution: The dilator(s) shall be replaced every 2 years or after 208 cleaning and disinfection cycles which correspond to twice a week for two years.

Always follow medical practice and national requirements regarding biohazard when disposing of a used medical device.

10. Reporting

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

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Atos and the Atos logo are trademarks of Coloplast A/S.
© 2025 Coloplast A/S. All rights reserved.
Atos Medical AB