PROVOX


Electrolarynx

TruTone Emote

TruTone Plus

SolaTone Plus

Oral adapter
Oral tube
Sound head

SolaTone
Plus

TruTone Emote

TruTone Plus

Figure 1


Figure 2


Figure 3


Figure 4



Figure 5

Limited Warranty

For international warranty, contact your local distributor.

Contents

EN - ENGLISH.....	6
DE - DEUTSCH.....	16
NL - NEDERLANDS.....	28
FR - FRANÇAIS.....	39
IT - ITALIANO	50
ES - ESPAÑOL.....	61
PT - PORTUGUÊS.....	72
SV - SVENSKA	83
DA - DANSK.....	93
NO - NORSK.....	103
FI - SUOMI.....	113
LV - LATVISKI.....	124
CS - ČESKY	134
HU - MAGYAR	144
SK - SLOVENČINA.....	155
PL - POLSKI.....	166
RO - ROMÂNĂ	178
HR - HRVATSKI.....	189
EL - ΕΛΛΗΝΙΚΑ.....	200
BG - БЪЛГАРСКИ.....	212
TR - TÜRKÇE	224
HY - ՀԱՅԵՐԵՆ.....	234
RU - РУССКИЙ.....	250
ID - BAHASA INDONESIA.....	263
MS - BAHASA MELAYU	274
JA - 日本語	285
KO - 한국어.....	298
ZHTW - 繁體中文.....	308
ZHCN - 简体中文	319
TH - ไทย.....	330
VI - TIẾNG VIỆT	343
HE - עברית	365
AR - العربية	375
FA - فارسی.....	386

1. Intended Use

An electrolarynx is a battery-powered artificial larynx that is applied externally and intended for use in the absence of the ability to use the anatomical larynx to produce sound.

When held against the skin in the area of the voice box, or by inserting the oral tube into the oral cavity (with an oral adapter), the device generates mechanical vibrations which resonate in the oral and nasal cavities and can be modulated by the tongue and lips in a normal manner, thereby allowing the production of speech.

Intended User Group

The patient is an intended operator.

2. Contraindications

The device should only be used in accordance with this IFU (Instructions for Use). Users without the physical, cognitive, or mental ability required to operate the device themselves should not use the device independently and should only use it if they are under sufficient supervision of a clinician or a trained caregiver. The device should not be directly applied over frail neck tissue with weak blood vessels, as this can cause tissue damage or bleeding. Patients with this condition should only use the device when they have been specifically instructed by their clinician on how to use the device and where to safely apply it.

3. Description of the Device

An electrolarynx is a battery-powered artificial larynx that is applied externally on undamaged skin and intended for use in the absence of the anatomical larynx or the inability to use the larynx to produce sound.

Device Types

TTE	Provox TruTone Emote
TTP	Provox TruTone Plus
TTE/P	Provox TruTone Emote/Plus
STP	Provox SolaTone Plus

Box Content

Electrolarynx	IFU (Instructions for Use)
Oral adapter	Oral tube variety pack
Micro-USB cable (1 m)	Extra sharp sound head (TTE/P)
AA NiMH rechargeable batteries x2 (STP)	Lanyard

Part Specification (see Figure 1)

A. Sound head (TTE/Px2) (STPx1) (including sound button and foam ring)	B. Sound head spacer (TTE/P)
C. Power button (Emotion™ button in TTE/P)	D. Volume thumbwheel
E. Body	F. Battery (STP)
G. Battery cap	H. Battery connector (TTE/P)
I. Oral tube (2+2)	J. Oral adapter
K. Pitch adjustment thumbwheel	L. LED indicator light
M. MODE button (TTE/P)	N. Battery (TTE/P)
U. USB dust cover (TTE/P)	V. USB charge jack

4. Warnings

- Incorrect use of the device or failure to follow this IFU can cause personal injury and/or damage to the device.
- Do not turn on, use, or charge the device if it has been exposed or dropped into water or other liquids, as this can cause personal injury and/or damage to the device.
- Use only approved chargers, as usage of unapproved chargers can cause fire, personal injury, or damage to the device.
- Before charging, inspect the charger for damage. A damaged or non-functioning charger shall no longer be used, as it can cause fire, personal injury, or damage to the device.
- Do not use the device while charging, as this can cause electric shock.
- Do not attempt to charge non-rechargeable batteries, as this could result in personal injury or damage to the device.
- Do not replace rechargeable batteries with non-rechargeable batteries, as this could result in personal injury or damage to the device.
- Do not physically damage the batteries, as this could result in personal injury or damage to the device.
- Do not store batteries in a pocket, purse, box, drawer, or similar where they may short-circuit each other, or be short-circuited by conductive materials such as coins or keys, as this could result in personal injury or damage to the device.
- Do not expose the lithium battery to temperatures outside the temperature ranges specified on the front of this IFU, as this could result in personal injury or damage to the device.
- Use caution when replacing the batteries. Incorrect installation of the batteries may result in a hazardous situation, such as extreme heat or fire, which could result in personal injury or damage to the device or other property.
- Do not use a leaking battery, as this could result in personal injury or damage to the device.
- Do not insert any part of the device or accessories into the stoma, as this could result in personal injury.
- Do not modify the device, as this could result in personal injury or damage to the device.
- Use of Electrolarynx near other electronic equipment should be avoided, as it could cause operational problems. If necessary, ensure that both devices are working normally.
- Portable RF (radio frequency) communications equipment, including peripherals such as antenna cables and external antennas, should be used no closer than 30 cm (12 in.) to the device, including cables specified by the manufacturer. Otherwise, it could affect the performance of the device.
- Use of cables other than those specified or provided with the device could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- There is a risk of strangulation due to long cable. Keep out of reach of children and pets.
- The device contains small parts that may become dislodged and pose a choking hazard. Keep away from small children.
- Do not use the device between persons. This can cause cross-contamination which can lead to serious infection. The device is intended for single patient use only.
- The device contains a magnet that generates magnetic and electromagnetic fields that may interfere with pacemakers or other implantable devices as well as certain procedures or treatments. Maintain a minimum distance of 15 cm (6 in.) between the device and any medically implanted devices. Consult with your physician before any medical procedure or treatment. If interference between the device and any medically implanted device is suspected, discontinue use and consult with your physician.

5. Precautions

- Be careful not to damage the battery wiring when opening/closing the battery cap.
- Failure to follow the maintenance instructions could result in personal injury and/or damage to the device.
- Do not charge the device from a PC as this may cause the battery to overheat and this could result in personal injury and/or damage to the device.
- Never use a broken, modified, or altered device; or charge a broken, modified, or altered device, as this could result in personal injury and/or damage to the device.
- The device and batteries may become warm, up to 45 °C (113 °F), while charging. Allow the device to cool down before use.
- The device may become warm during use. Take necessary precautions when placing it on your neck for use.
- The device may malfunction or stop working due to electromagnetic field interference from other devices.
- Be careful not to chew/bite on the oral tube, as this could cause dental damage.
- Be careful not to drop the device, as this could cause damage to the device.

Side-Effects

- With extensive use, the vibrations produced by the device could be perceived as uncomfortable by some users.
- The device may cause a sense of effort and fatigue in the shoulder/arm, especially in individuals with known shoulder problems.

6. How to Use

Note: Before use, allow the device to reach operating temperature (see Section 10 Technical Data). Let it cool or get warm for up to 2 h after transport or after storage at high or low temperatures.

Non-standardized Symbols (TTE/P)

Symbol	Symbol Label	Symbol Placement on the Device
	Pitch	Below the pitch adjustment thumbwheel
	Speaking	Below the power button
	Volume	Below the volume thumbwheel

6.1 Power On/Off using power button

- Power on by pressing the power button (C).
- Power off by releasing the power button.

Note: To avoid the device to be powered on and become hot during transportation in a bag or similar, turn the volume thumbwheel all the way down to turn off the volume before transport.

6.2 Adjusting the Volume and Power On/Off

1. Turn the volume thumbwheel (D)
 - all the way down to turn off volume and power.
 - up to turn on the device and adjust the volume.

The power button can also be used for pitch control by pressing the button firmly for a higher pitch and lightly for a lower pitch.

6.3 Adjusting the Pitch

1. Turn the pitch adjustment thumbwheel (K) to adjust the pitch.
2. Press and hold the power button (C) while making your first adjustment to hear the differences. Turn the thumbwheel up for higher pitch and down for lower pitch.

If you use the tone adjustment thumbwheel (K) without holding the power button (C), the tone will still change without you being able to hear it.

6.4 Setting the Mode (TTE/P)

1. To change the mode, first remove the battery cap and then press both the mode button (M) and the power button (C).

The mode changes immediately, but if you continue to hold down the power button, the device will beep to indicate the new mode setting. The number of beeps and the tone of the beep indicate which mode it is in.

Provox TruTone Emote (TTE) has six modes (1–6).

Provox TruTone Plus (TTP) has two modes (1–2).

Modes 1–4: Four preset ranges

Mode 1	Nearly monotone	1 beep
Mode 2	Low pitch variations (1/2-octave range)	2 beeps
Mode 3	Medium pitch variations	3 beeps
Mode 4	High pitch variations (2-octave range)	4 beeps

Modes 5–6 (TTE):

Volume mode allows you to control the volume with the power button (C) instead of the pitch. This is for people who do not use the pitch control, and it allows volume control without moving the thumbwheel. A very light touch lets you speak with your ‘whisper voice’, while pressing firmly emits sound at full volume. The volume thumbwheel still limits your maximum volume. Changing into these modes creates a ‘whoop’ sound (with an increasing tone), rather than a steady tone beep.

Mode 5	Low Sensitivity—full volume requires more pressure.	2 ‘whoops’
Mode 6	High Sensitivity—full volume requires less pressure.	3 ‘whoops’

6.5 Choosing Sound Head

If you have hard neck tissue or are in a loud environment, the grey sound head can be used because of its sharper sound (i.e. higher pitch) (see Section 7.2 Replacing the Sound Head).

6.6 Using the Sound Head Spacer (TTE/P)

If you regularly use a higher pitch and prefer the tonal qualities without a sound head spacer, you may remove it (see Section 7.3 Adding/Removing the Sound Head Spacer).

6.7 Charging the Battery (TTE/P)

1. Pull the USB dust cover (U) out at the top end.
2. Insert the micro-USB cable into the USB charge jack (V) using an approved wall charging plug* (not included).
3. Check that the LED indicator light (L) is lit.
4. Unplug the device when the charging is complete to save power. Use until volume begins to fade, which is 80% drained.

*approved Limited Power Source (LPS) according to IEC 60950-1 or IEC 62368-1

LED Indicator Light (TTE/P)

Low battery: The LED indicator Light (L) turns on when you press the power button, signaling that it is close to charging time.

Charging—Solid light
Charging completed—Flashing light

6.8 Placing Electrolarynx on the Neck

1. Place the sound head (A) against your neck and press the power button (C). Do not insert any part of the device or accessories into the stoma.
2. Make sure that the entire sound head is in contact with the skin of your neck, otherwise sound will escape and make a noise that makes it difficult to understand you.
3. Move your lips, mouth, jaw, and tongue (i.e. articulate) as you would normally

- speak. Press the button to start the sound, start articulating, and then release the button at the end of the phrase. Breath relaxed and natural—do not force the air out.
4. Try different positions until you find your "sweet spot" (i.e. the spot on your neck where the voice sounds the best). Note that even a 3 mm (1/8 in.) position change can have a great impact on the sound volume.

If you are unable to achieve sound transmission through your neck, or are unable to place the device against your neck for medical reasons, try cheek placement or use the oral adapter (see Section 7.4 Installing/Using Optional Oral Adapter/Tube). Adjusting the base pitch up or down may produce a better voice. Mid-range pitches are easier to hear for most people.

6.9 Optional Oral Adapter and Oral Tubes

If neck placement does not suit you due to a tender neck or excessive buzzing, you can use the oral adapter (see Section 7.4 Installing/Using Oral Adapter/Tube).

6.10 Guidelines for Demonstration Use (Valid for Health Care Professionals)

Warning: The oral adapter, oral tube, and sound head are for single patient use only and must not be reused between patients.

Precaution: Maintain the device after each patient and demonstration use according to the maintenance instructions.

- Demonstration use is intended to determine patient suitability for the device.
- The users shall clean and disinfect their hands with isopropyl alcohol (IPA) or with other suitable disinfectant, or wear gloves.
- The sound head, including the sound button and foam ring, shall be replaced and discarded between each patient (see Section 7.2 Replacing the Sound Head).
- The device shall be wiped off before and after each patient (see Section 7.1 Maintaining the Device).

7. Maintenance and Replacement of Part

7.1 Maintaining the Device

The user is responsible for the maintenance tasks. Perform maintenance after each use.

1. Wipe the device with a clean, dry cloth, or if necessary, a slightly damp (not dripping wet) cloth (see Figure 2). Use water with mild soap or isopropyl alcohol 50%–75%.

Precautions:

- Allow the device to dry before use.
- Be careful not to get moisture into the device.
- Do not use electrical or spray cleaners.

7.2 Replacing the Sound Head

Contact your local representative for replacement sound heads.

1. Wipe down the device (see Section 7.1 Maintaining the Device). Be sure to apply light pressure and wipe the power button with a circular or twisting motion. After wiping down the body, wipe down the sound head cap. Be careful not to apply excess pressure that could cause the top sound button to be dislodged.

2. Unscrew the sound head (A) and remove it (see Figure 2).

Note: Make sure not to press on the sound button when handling the sound head cap as this may cause the button and foam ring to be pushed through.

Make sure that the sound head spacer (B) is not broken or missing. Replace it, if necessary, before installing the new sound head (see Section 7.3 Adding/Removing the Sound Head Spacer).

3. Install the new sound head.

Caution: Do not remove the silicone diaphragm underneath the sound head. The silicone diaphragm should be retained in the groove on the coil/actuator stem (see Figure 3.1).

Caution: Do not twist the coil/actuator. The wires should wind smoothly around the coil/actuator (see Figure 3.2).

7.3 Adding/Removing the Sound Head Spacer

1. Unscrew the sound head (see Figure 4).
2. Add/remove the spacer.
3. Screw the sound head back on.
4. Test the device, listening to the tonal qualities. If you prefer the tonal qualities without the spacer, be sure to store the spacer in a safe place, such as in the box your device came in.

7.4 Installing/Using Oral Adapter/Tube

1. Insert the oral tube (I) into the top of the oral adapter (J).

Caution: Do not insert the oral tube beyond stop. Ensure that the tube does not protrude into the cap area (see Figure 3.1).

2. Place the oral adapter over the sound head (A).

Caution: Do not twist it on (see Figure 3.2).

3. Place the oral tube into the mouth, at the corner of the mouth.
4. Activate electrolarynx. Speak around the oral tube.

7.5 Changing the Battery

Valid for STP

Use two 1.5 V Nickel Metal Hydride (NiMH) rechargeable AA cells.

1. Install the batteries (F) by following the arrows placed inside the STP.

Valid for TTE/P

Use only approved Li-Ion batteries. Contact distributors for more information.

1. Install the battery (N). It only fits in one way. The battery connector (H) alignment tab points up, with the red wire on the right.

8. Troubleshooting Guide

8.1 No Satisfactory Intelligibility During Phone Use

- Place the phone microphone up by the nose, not below the mouth. This will get the microphone away from stoma and electrolarynx sound.
- Reduce the volume of the electrolarynx so that the volume level is as quiet as it can be, while still generating a tone.

8.2 Device is 'Dead', 'Broken', or 'Just stopped working'

- Make sure that the volume thumbwheel (D) is not turned all the way down (the device is turned off). Turn the volume thumbwheel (D) all the way down to off.
- Lightly press the power button (C) and increase pressure. If the device activates and deactivates in the process, the device needs to be sent in for replacement (see Section 9 Service and Assistance).
- If the LED turns on every time the power button (C) is pressed, the device is low on battery and needs to be charged (valid for TTE/P).
- Remove the battery cap and press the battery cable connector (H) further into the housing.

8.3 LED indicator Light Does not Turn on While Charging

- Test the power outlet with a lamp.
- Make sure the micro-USB cable is fully connected to the power supply and the device.
- Try a different micro-USB cable.
- Try a different power supply (USB power outlet).
- Test charger and cable on a cell phone or other device.

Valid for STP

- Make sure the batteries are rechargeable. The device does not charge non-rechargeable batteries.

- Make sure the batteries are inserted correctly. Check to make sure the batteries have proper orientation with the battery + going to the device +.

8.4 The device 'Has too much buzzing,' 'Sounds funny,' 'Is muffled,' or 'Not clear'

- Make sure the sound head (A) is in full contact with the skin on the neck.
- Try a different location on the neck (see Section 6.8 Placing Electrolarynx on the Neck).
- Reduce the volume (D).
- Adjust the pitch (K).
- Use the oral adapter (see Section 7.4 Installing/Using the Oral Adapter/Tube).
- Inspect the cap. If the sound button seems loose or has moved slightly, or if the foam ring seems damaged, the device might need a new foam ring.
- Remove the sound head (A) and look inside.
- Is there a small rubber piece glued in place in the center of the sound head (A)? If not, the device needs a replacement sound button (see Section 9 Service and Assistance).
- Make sure the diaphragm is attached to the coil stem. Make sure the coil bounces freely when touched, and does not have any gritty feel (see Figure 5).

Valid for TTE/P

Try the optional extra sharp sound head (see Section 7.2 Replacing the Sound Head).

8.5 'Sound head (A) came apart' or the 'Sound button has been pushed out of the cap'

Check if the foam ring is damaged. If it is not damaged, the original foam ring may be reused and reassembled. If it is damaged, it must be replaced.

Valid for TTE/P

The sensitivity of the button sensor may be too sensitive. Try a different MODE, usually a lower mode.

8.6 Button 'Sticks'

- Read the general cleaning instructions (see Section 7.1 Maintaining the Device).
- Use a Q-tip or cotton swab with a small amount of recommended alcohol to clean around the power button 5–6 times, then press the power button 5–6 times.
- Alternate between pressing and cleaning around the power button with a Q-tip or cotton swab. Continue until it no longer sticks.
- Use a thicker piece of paper to remove debris stuck in the notch around the power button.

8.7 Unsteady Tone (TTE/P)

Try a mode with a lower sensitivity (see Section 6.4 Setting the Mode).

8.8 Unsteady Volume (TTE/P)

Try a mode that is not a volume mode (see Section 6.4 Setting the Mode).

8.9 Device Causing Interference to Radio or TV Reception

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

9. Service and Assistance

The electrolarynx does not contain any serviceable parts other than those removable by the patient. Service parts include batteries, sound head, sound head spacer, foam ring, oral adapter, and oral tubes. Contact your local representative for assistance in using or maintaining your electrolarynx and to order available parts.

10. Technical Data

Description	Specification
Operating temperature (to maintain optimal battery life)	5 °C to 40 °C (41 °F to 104 °F) (valid for TTE/P) 5 °C to 25 °C (41 °F to 77 °F) (valid for STP)
Storage and transport temperature	-25 °C to 70 °C (-13 °F to 158 °F) (valid for TTE/P) -20 °C to 25 °C (-4 °F to 77 °F) (valid for STP)
Operating humidity (to maintain optimal battery life)	15%–90% relative humidity 700 hPa to 1060 hPa
Storage humidity	0%–45% relative humidity
Expected service life	3 years
Applied part	Type BF Applied Part, sound head
Mode of operation	TTE/P: Duty cycle 50%. For example, the device operates with an ON time of 1 min followed by an OFF time of 1 min. Please note that shorter ON times will require correspondingly shorter OFF times. STP: 30 min of activation during a 24-hour period
Power source	Internally powered
Dimensions	Provox SolaTone Plus 118 mm (4.6 in.) Provox TruTone Emote 113 mm (4.5 in.) Provox TruTone Plus 113 mm (4.5 in.)
Weight	Provox SolaTone Plus 130 g (0.29 lb) Provox TruTone Emote 130 g (0.29 lb) Provox TruTone Plus 131 g (0.29 lb)
IP Classification	IP22* (valid for TTE/P)
Charge specification	5 V, 750 mA minimum (or 0.75 A minimum)** max 100 VA

*The device was tested and remained safe during and after the test. However, the device may be temporarily inoperable or permanently damaged.

**Approved Limited Power Source (LPS) according to IEC 60950-1 or IEC 62368-1

Electromagnetic Emissions

Emission Test	Compliance	Electromagnetic Environment
RF emissions, CISPR 11	Group 1	The system uses RF energy only for its internal function. Therefore, its RF emission are low and are not likely to cause any interference in nearby electronic equipment.

Emission Test	Compliance	Electromagnetic Environment
RF emissions, CISPR 11	Class B	The system is suitable for use in all locations including those allocated for in residential environments and those directly connected to a low voltage power supply network which supplies buildings used for domestic purposes.
Harmonic emissions, IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emission, IEC 61000-3-3	Complies	

Note: The emission testing includes the use of an approved wall charging plug.

Electromagnetic Immunity (TTE/P)

Immunity Test	IEC 60601-1-2 Test Level	Compliance Level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±4 kV, ±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air (see also Electromagnetic Interference)
Electrical fast transient/burst IEC 61000-4-4	±2 kV a.c. power port 100 kHz repetition freq.	±2 kV a.c. power port 100 kHz repetition freq.
Surge IEC 61000-4-5	±0.5 kV, ±1 kV line-to-line	±0.5 kV, ±1 kV line-to-line

Immunity Test	IEC 60601-1-2 Test Level	Compliance Level
Voltage dips and interruptions IEC 61000-4-11	0% U_T , 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315°	0% U_T , 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315°
	0% U_T , 1 cycle at 0°	0% U_T , 1 cycle at 0°
	70% U_T , 25 cycles at 0°	70% U_T , 25 cycles at 0°
	0% U_T , 250 cycles at 0°	0% U_T , 250 cycles at 0°
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m
Proximity magnetic fields IEC 61000-4-39	8 A/m, CW at 30 kHz 65 A/m, DC 50% at 134.2 kHz, 7.5 A/m, DC 50% at 13.56 MHz	8 A/m, CW at 30 kHz 65 A/m, DC 50% at 134.2 kHz, 7.5 A/m, DC 50% at 13.56 MHz

Note: U_T is the a.c. mains voltage prior to application of the test level.

Immunity Test	IEC 60601-1-2 Test Level	Compliance Level
Conducted RF IEC 61000-4-6	3 Vrms and 6 Vrms in the ISM and amateur radio bands 150 kHz to 80 MHz 80% AM, 1kHz sine	3 Vrms and 6 Vrms in the ISM and amateur radio bands 150 kHz to 80 MHz 80% AM, 1kHz sine
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2,7 GHz 80% AM, 1kHz sine	10 V/m 80 MHz to 2,7 GHz 80% AM, 1kHz sine
Proximity fields from RF wireless communication port IEC 61000-4-3	27 V/m, PM 18 Hz at 385 MHz 28 V/m, PM 18 Hz at 450 MHz, 9 V/m, PM 217 Hz at 710 MHz, 745 MHz, 780 MHz 28 V/m, PM 18 Hz at 810 MHz, 70 MHz, 930 MHz 28 V/m, PM 217 Hz at 1720 MHz, 1845 MHz, 1970 MHz 28 V/m, PM 217 Hz at 2450 MHz 9 V/m, PM 217 Hz at 5240 MHz, 5500 MHz, 5785 MHz	27 V/m, PM 18 Hz at 385 MHz 28 V/m, PM 18 Hz at 450 MHz, 9 V/m, PM 217 Hz at 710 MHz, 745 MHz, 780 MHz 28 V/m, PM 18 Hz at 810 MHz, 70 MHz, 930 MHz 28 V/m, PM 217 Hz at 1720 MHz, 1845 MHz, 1970 MHz 28 V/m, PM 217 Hz at 2450 MHz 9 V/m, PM 217 Hz at 5240 MHz, 5500 MHz, 5785 MHz

Note: The immunity testing includes the use of an approved wall charging plug.

Electromagnetic Interference

While exposed to electromagnetic discharge (ESD), the battery might need to be removed and put back again to restore normal operation.

11. Disposal

Always follow medical practice and national requirements regarding biohazards when disposing of a used medical device.

12. Reporting of Incidents

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

13. Compliance

The device is tested and conforms to the IEC 60601-1-2 with test level adopted for domestic environment (valid for TTE/P).

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UK Responsible Person
Atos Medical UK Ltd Tottle Road
Cartwright House
Nottingham
Nottinghamshire NG2 1RT
England United Kingdom

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Atos

atosmedical.com

Australian Sponsor:
Coloplast Pty Ltd
Level 7, 1 Peters Avenue,
Mulgrave, Victoria 3170
Australia
www.coloplast.com.au
+61 1800 286 728

Brazil Registration Holder:
Atos Medical Brasil Ltda
Rua Joaquim Nabuco, 47
Cj. 92- Brooklin
04621-000 Sao Paulo, SP
Brasil
Ph: +55 11 2305 2022

Canadian Distributor:
Coloplast Canada Corp.
Canadian Distribution Centre
1380 Creditstone Rd
Unit 6&7
Vaughan, ON L4K 0J1
Canada
+1 905 829 7651

Atos Medical AB
Kraftgatan 8, SE-242 35 Hörby, Sweden
Tel: +46 (0)415 198 00 • info@atosmedical.com

Atos and the Atos logo are trademarks of Coloplast A/S.
© 2025 Coloplast A/S. All rights reserved.