


Tracoe Twist Plus to Shiley Flex Comparison Chart

Tracoe Twist Plus tubes have a longer tube length than Twist tubes and an optimized inner to outer diameter ratio

Comparison Chart

The chart below is intended as a reference when evaluating conversion from a Shiley Flex to a Tracoe Twist Plus tracheostomy tube. Tracheostomy tube sizing and dimensions vary by manufacturer, please refer to individual product specifications when making tube selections.



Shiley Flex

Size	ID-IC	OD-OC	TL
Shiley Flex 4	5.5	9.4	62
Shiley Flex 5	6.0	10.1	68
Shiley Flex 6	6.5	10.8	74
Shiley Flex 7	7.0	11.4	77
Shiley Flex 8	7.5	12.2	79
Shiley Flex 9	8.0	12.7	79
Shiley Flex 10	9.0	13.8	79

Tracoe Twist Plus

Size	ID-IC	OD-OC	TL
Tracoe Twist Plus 7	7.0	9.8	85
Tracoe Twist Plus 8	8.0	10.8	88
Tracoe Twist Plus 9	9.0	11.8	90
Tracoe Twist Plus 10	10.0	12.8	92

All sizes are in mm ; **ID-IC:** Inner diameter of inner cannula; **OD-OC:** Outer diameter of outer cannula; **TL:** The length measured along the center line of the tube, starting behind the neck flange to the distal end

Ordering information

Size	Cuffed	Uncuffed
Tracoe Twist Plus 7	REF 311-07	REF 313-07
Tracoe Twist Plus 8	REF 311-08	REF 313-08
Tracoe Twist Plus 9	REF 311-09	REF 313-09
Tracoe Twist Plus 10	REF 311-10	REF 313-10


In sterile package and includes: inner cannulas with 15mm connector, perforated obturator, adjustable neck strap and product info card with three detachable labels. To order, specify: REF + Size, e.g. REF 311-07. Tracoe Twist Plus is also available in Fenestrated and Extract

Tracoe Twist to Shiley Flex Comparison Chart

Tracoe Twist tubes have a shorter tube length and smaller size options than Twist Plus tubes

Comparison Chart

The chart below is intended as a reference when evaluating conversion from a Shiley to a Tracoe Twist tracheostomy tube. Tracheostomy tube sizing and dimensions vary by manufacturer, please refer to individual product specifications when making tube selections.



Shiley Flex

Size	ID-IC	OD-OC	TL
Shiley Flex 4	5.5	9.4	62
Shiley Flex 5	6.0	10.1	68
Shiley Flex 6	6.5	10.8	74
Shiley Flex 7	7.0	11.4	77
Shiley Flex 8	7.5	12.2	79
Shiley Flex 9	8.0	12.7	79
Shiley Flex 10	9.0	13.8	79

Tracoe Twist

Size	ID-IC	OD-OC	TL
Tracoe Twist 4	4.0	7.2	59
Tracoe Twist 5	5.0	8.6	66
Tracoe Twist 6	6.0	9.2	72
Tracoe Twist 7	7.0	10.4	74
Tracoe Twist 8	8.0	11.4	76
Tracoe Twist 9	9.0	12.5	78
Tracoe Twist 10	10.0	13.8	80

Ordering information

Size	Cuffed	Uncuffed
Tracoe Twist 4	REF 301-04	REF 303-04
Tracoe Twist 5	REF 301-05	REF 303-05
Tracoe Twist 6	REF 301-06	REF 303-06
Tracoe Twist 7	REF 301-07	REF 303-07
Tracoe Twist 8	REF 301-08	REF 303-08
Tracoe Twist 9	REF 301-09	REF 303-09
Tracoe Twist 10	REF 301-10	REF 303-10

In sterile package and includes: inner cannulas with 15mm connector, perforated obturator, adjustable neck strap and product info card with three detachable labels. To order, specify: REF + Size, e.g. REF 301-04. Tracoe Twist is also available in Fenestrated and Extract

Atos Medical Inc

5000 S. Towne Drive, Suite 200, New Berlin, WI 53151-7956 USA
+1 800 217 0025 | info.us@atosmedical.com | www.atosmedical.us

Tracoe Twist

