



Tracoe Medical GmbH | Reichelsheimer Str. 1 / 3 | 55268 Nieder-Olm

TRACOE medical GmbH

Reichelsheimer Straße 1 / 3
55268 Nieder-Olm

Tel: 06136 9169-0

Fax: 06136 9169-200

info@tracoe.com

www.tracoe.com

Geschäftsführer:

Dr. Hans Kränzler

Amtsgericht Mainz, HRB 44548

Steuer-Nr.: 08/658/71878

Ihr Zeichen

Unser Zeichen

Datum

14. März 2024

To whom it may concern

MRI Safety Information

We, TRACOE medical GmbH, hereby declare that the MRI Safety Information of our Tracheostomy Tubes listed below is as follows:

REF	MRI Safety Information
Tracoe Vario, Tracoe Vario XL	
REF 450	Unsafe
REF 455	Unsafe
REF 460	Conditional
REF 464	Safe
REF 470	Conditional
REF 451	Unsafe
REF 461	Conditional
REF 471	Conditional
Tracoe Twist Plus	
REF 311	Conditional
REF 312	Conditional
REF 313	Safe
REF 314	Safe
REF 316	Conditional
REF 888-316	Conditional

Commerzbank Hanau
IBAN DE63 5064 0015 0240 0604 00
BIC COBADEFF

Frankfurter Volksbank
IBAN DE48 5019 0000 6100 0113 00
BIC FFVBDEFF

Deutsche Bank Frankfurt
IBAN DE83 5007 0010 0469 4006 00
BIC DEUTDEFF

Tracoe Twist	
REF 301	Conditional
REF 302	Conditional
REF 303	Safe
REF 304	Safe
REF 306	Conditional
REF 888-306	Conditional
Tracoe Mini	
REF 350	Safe
REF 355	Safe
REF 356	Unsafe
Tracoe Silcosoft	
REF 360	Conditional
REF 361	Conditional
REF 362	Conditional
REF 363	Conditional
REF 370	Conditional
REF 371	Conditional
REF 372	Conditional
REF 373	Conditional
Tracoe Comfort, Tracoe Comfort XL, Tracoe Comfort Plus	
Outer cannula	Safe
Inner cannula with low profile	Safe
Inner cannula with 15 mm connector	Safe
Inner cannula with 22 mm connector	Safe
Inner cannula with speaking valve type A/B	Unsafe

Yours sincerely,



Felix Dreher

Head of Product Lifecycle Management

Mitglied im

