

Tracoe Vario XL Extract – REF 471-P

Tracheostomy tube with minimally traumatic insertion system, adjustable neck flange, low-pressure cuff, subglottic suction, scale and 15 mm connector

Product Description

Tracoe Vario is a soft, flexible tube with an adjustable neck flange. By means of a practical push-button mechanism, the flange can be moved and individually adjusted for each patient. In this way, the correct length can be achieved. A scale ensures the exact positioning of the neck flange in preparation for a tube change. The flange is then locked in place by moving the orange lever upward. Locking with a twist lock is no longer necessary. Tracoe Vario tubes are made of DEHP-free, medical grade PVC and are available with and without spiral-reinforcement. The assortment includes both standard and extra-long versions.

Scope of Supply

- Obturator
- Adjustable neck strap
- Product info card with two detachable labels



General Information

Product line	Vario XL Extract
Legal manufacturer	Tracoe Medical GmbH Reichelsheimer Str. 1/3 55268 Nieder-Olm, Germany
Country of origin	Germany
Colour Coding	REF 471-P XL
Handling	  
Shelf life	5 years, from date of production
Maximum Useful Life	29 days, from first use
Maximum Useful Life Insertion System	single-use disposable product

Packaging

Packaging unit	1 piece
Weight	See Table
Dimensions (W x H x L)	271 x 47 x 113 mm
Material outer packaging:	Cardboard box
Material sterile packaging:	Sterile bag with medical paper and transparent heatsealable PET/PP foil

Material Specifications

Material Outer cannula	Polyvinylchloride
Material Neck Flange	Polysulfone/ TPE
Material Cuff	Polyvinylchloride
Biocompatibility	EN ISO 10993-1
   	

Storage




Certification

Product classification acc. MDR	IIb
CE Marking	 0 1 2 3
GMDN	35404
Quality system certified to	EN ISO 13485
FDA registration	Clearance Number K203362

Product Range									
Size	Order Number	ID	OD	TL (TL Min-TL Max)	OB	BA	CD	Gross weight	Global Trade Item Number
		mm	mm	mm	mm	°	mm	g	
08	REF 471-08-P	8	11.2	73.0 – 113.0	128	100	28	235	04035324019335
09	REF 471-09-P	9	12.3	73.5 – 118.5	134	100	30	236	04035324019342

Legend

- ID Inner diameter at bottom of cannula
- OD Outer diameter at bottom of cannula
- TL MIN Minimum length along centre line
- TL MAX Maximum length along centre line
- OB Length along outer bend
- BA Bending angle
- CD Cuff diameter



Tracoe Medical GmbH
Reichelsheimer Str. 1 / 3
55268 Nieder-Olm
Germany



Breathing-Speaking-Living
Atos Medical AB
Hyllie Boulevard 17
215 32 Malmö
Sweden
www.atosmedical.com