

Tracoe Twist – REF 302

Tracheostomy tube with low pressure cuff and fenestration

Product Description

The Tracoe Twist tracheostomy tubes are used as ventilation and therapy tubes with inner cannulas in hospitals, rehabilitation, and community settings. A special feature is the curved neck flange that moves in both the vertical and horizontal axes, following the movements of the patient's head. The thin tube walls of inner and outer cannulas allow for a good inside - outside diameter ratio. Tracoe Twist and Twist Plus tracheostomy tubes are made of polyurethane (PU). The low-pressure cuff is made of medical grade PVC.

Scope of Supply

- inner cannula with 15 mm connector
- inner cannula fenestrated with 15 mm connector
- Perforated obturator
- Neck strap



General Information	
Product line	Twist
Legal manufacturer	Tracoe Medical GmbH Reichelsheimer Str. 1/3 55268 Nieder-Olm, Germany
Country of origin	Germany
Colour Coding	REF 302
Handling	 
Shelf life	5 years, from date of production
Maximum Useful Life	29 days, from first use

Packaging	
Packaging unit	1 unit
Weight	See Table
Dimensions (W x H x L)	113 x 210 x 50 mm
Material outer packaging:	Cardboard box
Material sterile packaging:	Sterile bag with medical paper and transparent heatsealable PET/PP foil

Material Specifications	
Material Outer cannula	Polyurethane
Material Neck Flange	Polycarbonate
Material Cuff	Polyvinylchloride
Biocompatibility	EN ISO 10993-1
   	

Storage
  


Certification	
Product classification acc. MDR	IIb
CE Marking	
GMDN	35404
Quality system certified to	EN ISO 13485
FDA registration	Clearance Number K961449

Product Range

Size	Order Number	ID-IC	OD-OC bc	OD-OC bf	TL	OB	BA	CD	Gross Weight	Global Trade Item Number
		mm	mm	mm	mm	mm	°	mm	g	
04	REF 302-04	4,0	7,2	8,6	59,0	63,0	110	15,0	113	4035324033317
05	REF 302-05	5,0	8,6	10,1	66,0	73,0	100	18,0	115	4035324033416
06	REF 302-06	6,0	9,2	10,8	72,0	79,5	95	23,0	115	4035324033515
07	REF 302-07	7,0	10,4	12,0	74,0	80,5	95	26,0	116	4035324033614
08	REF 302-08	8,0	11,4	12,7	76,0	83,0	90	28,0	117	4035324033713
09	REF 302-09	9,0	12,5	14,2	78,0	88,5	90	30,0	123	4035324063314
10	REF 302-10	10,0	13,8	15,2	80,0	89,5	90	32,0	125	4035324033812

Legend

ID-IC	Inside diameter at bottom of inner cannula
OD-OC bc	Outside diameter at bottom of outer cannula
OD-OC bf	Outside diameter of outer cannula behind the flange
TL	Length along centre line from start of neck flange to bottom of tube
OB	Length along outer bend from start of neck flange to bottom of tube
BA	Bending angle
CD	Cuff diameter



Tracoe Medical GmbH
Reichelsheimer Str. 1 / 3
55268 Nieder-Olm
Germany



Breathing-Speaking-Living

Atos Medical AB
Hyllie Boulevard 17
215 32 Malmö
Sweden
www.atosmedical.com