

Tracoe Twist Plus – REF 521-XXX

Tracoe Twist Plus spare inner cannulas

Product Description

The Tracoe Twist Plus tracheostomy tubes are used as ventilation and therapy tubes with inner cannulas in hospitals, rehabilitation, and community settings. A special feature is the curved neck flange that moves in both the vertical and horizontal axes, following the movements of the patient's head. The thin tube walls of inner and outer cannulas allow for a good inside - outside diameter ratio. Tracoe Twist and Twist Plus tracheostomy tubes are made of polyurethane (PU). The low-pressure cuff is made of medical grade PVC.



Scope of Supply

- 30 inner cannulas with 15 mm connector

General Information	
Product line	Twist Plus
Legal manufacturer	Tracoe Medical GmbH Reichelsheimer Str. 1/3 55268 Nieder-Olm, Germany
Country of origin	Germany
Colour Coding	n. A.
Handling	
Shelf life	5 years, from date of production
Maximum Useful Life	5 weeks, from first use

Material Specifications	
Material Outer cannula	Polyurethane
Material Neck Flange	n. A.
Material Cuff	n. A.
Biocompatibility	EN ISO 10993-1
	

Storage


Packaging	
Packaging unit	30 units
Weight	See Table
Dimensions (W x H x L)	320 x 440 x 0,07
Material outer packaging:	Pressure lock bag
Material sterile packaging:	Sterile bag with medical paper and transparent heatsealable PET/PP

Certification	
Product classification acc. MDR	IIb
CE Marking	
GMDN	35404
Quality system certified to	EN ISO 13485
FDA registration	Clearance Number K961449

Product Range				
Size	Order Number	ID-IC	BA	Global Trade Item Number
		mm	°	
07	REF 521-07-XXX	7,0	100	04035324027989
08	REF 521-08-XXX	8,0	100	04035324027996
09	REF 521-09-XXX	9,0	100	04035324028009

Legend

- ID-ICInside diameter at bottom of inner cannula
- BA Bending angle



Tracoe Medical GmbH
Reichelsheimer Str. 1 / 3
55268 Nieder-Olm
Germany



Atos Medical AB
Hyllie Boulevard 17
215 32 Malmö
Sweden
www.atosmedical.com