


Tracoe Percutaneous Dilatation – REF 517

Tracoe Seldinger Guidewire - with guiding catheter

Product Description

Percutaneous Dilation Tracheostomy (PDT) has become routine in intensive care. The Ciaglia technique is now the method most widely used. The Tracoe Experc system is based on this proven technique. The products are tailored to the use of fenestrated and non-fenestrated tracheostomy tubes, including tubes with subglottic suction. They are supplied already mounted on their minimally traumatic insertion systems (P-tubes).



General Information

Product line	Percutaneous Dilatation
Legal manufacturer	Tracoe Medical GmbH Reichelsheimer Str. 1/3 55268 Nieder-Olm
Country of origin	Germany
Handling	 
Shelf life	5 years, as from date of production

Packaging

Packaging unit	1 Unit
Weight	74 g
Dimensions (W x H x L)	113 x 272 x 61 mm
Material outer packaging:	Cardboard box
Material sterile packaging:	Sterile bag with medical paper and transparent heatsealable PET/ PP foil

Material Specifications

Material Guid Wire	Nitinol, TFE coated
Material Guiding Catheter	Flurinated ethylene propylene
Biocompatibility	EN ISO 10993-1



Storage



Certification

Product classification acc. MDR	IIa
CE Marking	 0123
GMDN	11263
Quality system certified to	EN ISO 13485
FDA registration	Clearance number K043160 / D020866
Global Trade Item Number	04035324016310



Tracoe Medical GmbH,
Reichelsheimer Str. 1 / 3,
55268 Nieder-Olm
Germany



Breathing-Speaking-Living

Atos Medical AB
Hyllie Boulevard 17
215 32 Malmö
Sweden
www.atosmedical.com